

Reopen with confidence

Focus on reopening safely with an easy-to-implement COVID-19 testing experience

The BD Veritor™ At-Home COVID-19* Test enables reliable, self-administered COVID-19 testing with digitally displayed, shareable results in just 15 minutes.

*EUA Authorized by FDA

Veritor™
RAPID DIGITAL TEST
&Go

BD Veritor™
At-Home COVID-19 Test



Rapid, unambiguous results

- Clear, digitally displayed results in 15 minutes
- Trusted hospital-grade technology
- Results delivered via the Scanwell® Health App



Results you can trust

- Analyzes, interprets and displays clear, easy-to-read digital results in words, not lines, which minimizes the risk of incorrect result interpretation
- Internal control mechanism to confirm, without the presence of a healthcare provider, that a human sample is collected
- Results are verified by the app, not interpreted and entered by the user, increasing confidence in results



Ease of use

- Pain-free nasal sample collection
- Step-by-step video instructions and timers provided in the Scanwell® Health App
- Authorized for non-prescription, home use with self-collected (unobserved) nasal swab specimens from individuals aged 14 and older




Seamless result sharing

- Stores a secure digital record of results that can be accessed at any time
- Automated reporting to public health authorities
- Ability to integrate with digital platforms that manage employee health status at scale through customized dashboards, reports and alerts

Confidence in at-home COVID-19 testing

The **BD Veritor™ At-Home COVID-19* Test** offers seamless integration for rapid COVID-19 testing from the comfort of your home.

Simple in-app workflow guides user through the testing process




SIGN UP/LOG IN

SELECT PROFILE

STEP-BY-STEP INSTRUCTIONS

SCAN




RESULTS SCREEN



Automatic reporting to health department




If your organization has opted into result sharing, users can securely share results through the app



OPTIONAL RESULT SHARING

BD is constantly working to update compatibility and add more smartphones. Check if your smartphone is compatible and get the Scanwell® Health App



BDVeritorAtHome.com/devices

The **BD Veritor™ At-Home COVID-19 Test** What's in the box?

1. 1 kit box with tube holder
2. 2 nasal swabs
3. 2 scan cards
4. 2 test sticks
5. 2 tubes (with liquid) in a foil pouch
6. Quick Start Guide
7. Product Information Leaflet
8. Fact Sheet for Individuals



*This product has not been FDA cleared or approved; but has been authorized by FDA under EUA. This product has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens

The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of IVDs for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

Click or scan for
more information
or to order



BDVeritorAtHome.com

BD Life Sciences, 7 Loveton Circle, Sparks, MD 21152-0999 USA
800-638-8663

BD, the BD Logo, and Veritor are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2021 BD. All rights reserved. (BD-46120 1318-US-1121 November 2021)

